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EN | INSTRUCTIONS FOR USE

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ADITAMENTOS, COMPONENTES
E INSTRUMENTAL

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
EN | INSTRUCTIONS FOR USE

Abutments, components and instruments

INSTRUCTIONS FOR USE

IMPORTANT INFORMATION

Read these instructions for use carefully before using the product.

This document contains basic information for the use of the prosthetic instruments and abutments or components of the ZIACOM® original dental implant systems, from now on referred to as “ZIACOM® products”. This document has been written as a quick consultation guide for the practitioner responsible for the treatment, from now on referred to as “User”, and does not provide sufficient technical specifications for the correct use of the ZIACOM® products. It is not an alternative or substitute for specialised training and professional clinical experience.

The ZIACOM® products must be used carrying out sufficient planning of the treatment and rigorously following the surgical and/or prosthetic protocols set by Ziacom Medical SLU. Read the specific surgical and prosthetic protocols for each product carefully, as well as the instructions for use and maintenance, before using any ZIACOM® product. You can read them on our website, www.ziacom.es, or request them from your official authorised Ziacom Medical SLU distributor.

1. Description of the system

The ZIACOM® products are composed of dental implants, prosthetic components or components and surgical and prosthetic instruments.

Ziacom Medical SLU has created a wide range of dental implant systems for the carrying out of dentistry treatments in accordance with the existing scientific bibliography and current clinical standards. Ziacom Medical SLU develops, manufactures and commercialises latest-generation implant products and a wide range of restorative solutions to respond to the particular needs of each patient.

2. Description of the products

Instruments

The ZIACOM® instruments are manufactured in stainless steel, grade 5 titanium and/or plastic. Some instruments are laser-engraved to facilitate identification and use.

ZIACOM® instruments

Surgical drills – Surgical taps – Insertion keys – Circular scalpels – Extenders – Surgical Trephine Drills – abutment adaptors – Screwdrivers – Drill stops – Bone expanders – Osteotomes – Ratchet wrenches – Screwdriver tips – Thread adaptors – Retreatment positioners – Paralleling pins – Zinic®3D rings – Depth measurers – Abutment extraction screws – Radiographic templates.

Abutments

The ZIACOM® prosthetic components or components are manufactured in grade 5 titanium, Cr-Co alloy and biocompatible plastic materials. Some ZIACOM® components, abutments and screws have an anodised colour depending on the type of prosthetic platform to facilitate identification and use (see catalogue).

ZIACOM® prosthetic components and components

Preparable straight abutments – Preparable angled abutments – Transepithelial abutments – Abutments for overdenture – Interphase abutments – Provisional abutments – Clinical screws – Cover screws – Healing abutments – Abutments with machined base – Impression abutments – Laboratory test abutments – Impression screws – Laboratory screws – Analogues – Scanbody – UCLA.

3. Indications for use

It is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and study of models, as well as to determine the possible contraindications for the treatment, prior to carrying out a treatment with ZIACOM® products. A correct treatment with dental implants must be based on a diagnosis, a treatment plan and a design of the prosthetic rehabilitation which meets the needs and expectations of the patient.

Instruments

The ZIACOM® surgical and prosthetic instruments are designed as one of the fundamental parts of the implant treatment. With the instruments, you will be able to carry out all the necessary processes for the preparation of the bone tissue, the insertion of the ZIACOM® dental implants and the manipulation of the different ZIACOM® prosthetic components or components.

The surgical drills and trephine drills are designed to withstand a maximum of 45 uses under optimum conditions, whilst the circular scalpels are designed to withstand a maximum of 25. Improper use, lack of cleaning and maintenance or inadequate disinfection can reduce its number of uses.

The ZIACOM® instruments are re-usable.

Abutments

The ZIACOM® prosthetic components or components are fixed to the ZIACOM® dental implants and serve as a support for the prosthetic rehabilitation. The prosthetic rehabilitation may be single or multiple, partial or total, and may be screwed or cemented to the abutments, or may be removable by means of implant anchorages or retainers.

The attachments or abutments in biocompatible plastic material (PEEK) are for provisional use and are designed to be used in the mouth for a period of time no longer than 180 days. The removable prostheses with anchorages or retainers must be monitored at least once per year by the clinical professional. The correct intraoral adjustment of the prosthesis and the retention between the anchorages or retainers and the dental prosthesis must be verified. A loss of adjustment and/or a poor retention of the prosthesis put the prosthetic rehabilitation and the dental implants in danger.

The ZIACOM® abutments are single-use products.

4. Contraindications

It is necessary to carry out a preoperative medical examination of the patient and to determine if absolute and/or relative contraindications or risk factors for the treatment with dental implants exist. The user is responsible for evaluating the potential risks and benefits of the implant treatment and for taking the pertinent decisions depending on the case. The ZIACOM® products must not be used on patients who lack the minimum medical conditions for receiving a prosthetic treatment and rehabilitation with dental implants.

Relative contraindications

Age of the patient- Uncontrolled systemic disease- Poor oral hygiene- Active local infection - Pregnancy- Therapy with bisphosphonates- Sjogren syndrome- Uncontrolled endocrine disease - Radiation to head and neck- Allergy to stainless steel.- Stress - Psychological factors- Valve prostheses- Uncontrolled periodontal disease- Smoking- Alcoholism- Drug addiction- Addictions in general.

Absolute contraindications

The manufacturer does NOT recommend the implant treatment in the following cases: Cardiovascular diseases or disorders – Haemophilia – Cirrhosis of the liver – Terminal pathologies – Epilepsy – Neoplasms under care with chemotherapy – Blood dyscrasias – Cancer – Uncontrolled diabetes – Untreated periodontitis – Infectious diseases – Leukaemia – Alterations of the immune system – Treatment with immunosuppressives – Allergy to titanium – Active treatments with radiotherapy – Indefinite treatments with corticosteroids – Treatments with anti-coagulants – Involuntary muscle disorders – Serious liver disorders – Endocrine disorders – Psychiatric and/or psychological disorders.

Ziacom Medical SLU will update these contraindications when necessary, based on the information obtained through clinical studies and post-sales information.

5. Storage

The ZIACOM® products must be stored in their original packaging, in a clean, dry place, protected from sunlight and adverse conditions.

After cleaning and sterilisation, the torque wrench must be stored in the relaxed position at approximately 10 Ncm.

6. Packaging and supply state

The ZIACOM® instruments and prosthetic components or components are supplied packaged in a thermosealed plastic bag; the product inside is identified with a label stuck to the back.

The torque wrench is supplied in the 10 Ncm position and with an accuracy of +/- 10% of the set torque value.

Ziacom Medical SLU recommends that the traceability of the abutments used be maintained by recording the reference and the lot used in in the patient’s file.

Ziacom Medical SLU guarantees that all its products follow an exhaustive process of manufacture, control and cleaning before being packaged. The correct integrity of the packaging must be controlled before using the ZIACOM® products, checking that it is not damaged. If there is any defect or damage, it must not be used and must be notified to Ziacom Medical SLU or to the official authorised distributor.

7. Cleaning and sterilisation

The instruments and prosthetic components or components of the ZIACOM® dental implants systems are supplied NON-sterile.

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[For USA, the recommend validated sterilization cycle using steam and gravity autoclave is: 132°C / 270°F, minimum 15 minutes and 15-30 minutes Dry Time].

8. Precautions

Surgery and Instruments

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the dental implants. Prior to the intervention, it is essential to have consulted the surgical protocol of the product, as well as all the indications about its correct usage indicated in the catalogues of the ZIACOM® products.

Particular attention must be paid and the necessary measures taken to prevent the instruments, of such small size, from being swallowed or breathed into an airway. All the ZIACOM® instruments for manual use incorporate retaining elements between each other or for use with dental floss, to avoid accidental ingestion.

It is important to reduce the time of surgery and the damage done to the area receiving the implants, paying special attention to bone trauma and the asepsia of the surgical field. The correct preparation of the surgical site is a fundamental factor in the success of the treatment. This process requires that the surgical protocol for drilling set by the manufacturer for each type of implant and the correct use of the surgical drills, both in terms of the number of uses and of the working speeds, be followed rigorously. The surgical drills and taps must be used with constant abundant irrigation to prevent overheating of the bone.

Take special care not to apply excessive forces or torques during the insertion of the implant. It is important to control the insertion torque of the implant, both in insertion with contra-angle, and with insertion with ratchet wrench.

Attention

When a surgical instrument with a limited number of uses is used, that number must not be passed, as this would put the process of osseointegration, and therefore the treatment, in grave danger.

The maximum torque for insertion of the dental implants is 50 Ncm. Exceeding the maximum insertion torque indicated for the implants can cause serious damage to the dental implant, to its connection, to the transporter and to the clinical screw included.

Prosthetic rehabilitation

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the prosthetic rehabilitation. The design, the type of rehabilitation and the execution of the dental prosthesis must be studied prior to the insertion of the dental implants. Take special care not to apply excessive forces during the tests of the abutments or dental prostheses; this must always be done manually. Pay special attention to the definitive adjustment of the rehabilitation; the torque applied to the definitive clinical screws must be correctly controlled following the manufacturer’s indications, whether carried out mechanically or manually.

ZIACOM® Torque table for abutments:

Element/Abutment	Instrument/Tool	Torque
Cover screws/healing abutments	Hex. screwdriver 1,25mm	Manual
Impression abutment screws	Hex. screwdriver 1,25mm	Manual
Laboratory screws	Hex. screwdriver 1,25mm	10 Ncm
Clinical screws direct to implant	Hex. screwdriver 1,25mm	30 Ncm
Kiran® clinical screws direct to implant	Hex. screwdriver 1,25mm	30 Ncm
Basic®/XDrive® abutments	Insertion key: MABA100/MABA200	30 Ncm
Clinical screws on Basic®	Hex. screwdriver 1,25mm	25 Ncm
Kiran® clinical screws on Basic®	Hex. screwdriver 1,25mm	25 Ncm
Clinical screws on XDrive®	Hex. screwdriver 1,25mm	20 Ncm
Kiran® clinical screws on XDrive®	Hex. screwdriver 1,25mm	20 Ncm
ZiaCam® Scanbody + screw	Hex. screwdriver 1,25mm	Manual
Kirator® abutments	Insertion key: LOSD01/LOSD02	30 Ncm
Equator® abutments	Hex. screwdriver 1,25mm	30 Ncm
Tx30® abutment/Screw (Variable Rotation)	Torx.Tx30®screwdriver	30 Ncm

Some instruments are designed for use with a contra-angle. In order to guarantee a correct connection of both products, please see the contra-angle manufacturer’s instructions for use.

Attention

Applying insufficient torque can cause incidents of loosening of said screws. If this occurs, we recommend replacing them with new products, as retightening a screw increases the risk of fatigue and fracture of it.

Exceeding the maximum torque applied to the clinical screws indicated by the manufacturer can cause serious damage to the dental implant, to its connection, to the prosthetic abutment and to the clinical screw. Always check the passive adjustment of the implant-prosthesis unit. An incorrect passive adjustment puts the prosthetic rehabilitation in grave danger and can cause future damage to the dental implant, to its connection, to the prosthetic abutment, to the clinical screw and/or cause the loss or reabsorption of bone.

9. Warnings

The planning of the treatment and the fitting of the ZIACOM® dental implants requires specific dentistry training. The specifications of the product alone do not guarantee a proper use of it. The users are recommended to attend specialised theoretical/practical training courses to become familiar with the work techniques and protocols suitable for each product, including biomechanical, radiographical and prosthetic requirements associated with the implant treatment. It is the user’s obligation to inform themselves of and familiarise themselves with the current state of technique for each product and of its possible applications.

Before using ZIACOM® products, it is necessary to know and be familiar with the corresponding surgical and prosthetic procedures. The patient must meet sufficient anatomical and psychological conditions to undergo the treatment with dental implants. The user must ensure the correct planning of the treatment, guaranteeing its safety with a minimum margin for error and respecting the vital oral structures and the general health of the patient.

The procedure for use of the surgical drills, taps and other surgical instruments needed for the fitting of the implants is specified in the surgical protocols corresponding to each type of implant. The fitting of the implant and the prosthetic planning must adapt to the patient’s individual conditions, and in particular to a correct distribution of forces. The passive adjustment in the prosthetic rehabilitation structure and the adjustment of the occlusion must be achieved, and the appearance of excessive lateral forces avoided. A treatment with an insufficient number of implants, or an unsuitable choice of size and position to withstand and transmit the loads expected, can cause mechanical failure of the implant, of the abutment or prosthetic components or of the clinical screws due to overburdening or fatigue and substantial bone loss around the implant.

A qualitatively or quantitatively insufficient alveolar bone, the appearance of infections or diseases in general, and changes in the habits of the patient, are some potential causes of failure of the osseointegration and of the treatment. The lack of bone tissue or soft tissue can give rise to a deficient aesthetic result. An unsuitable prosthetic rehabilitation can cause the failure of the whole rehabilitation treatment.

The reuse of single-use products carries the danger of possible deterioration of their properties. The geometry of the product may have suffered damage during previous use or due to an unsuitable use. This implies the risk of surgical or prosthetic failure and/or harm to the patient’s health.

Electrosurgery is not indicated due to the conductivity of the dental implants.

10. Adverse or side effects

Many adverse or undesired effects can present themselves in current implantology. These are documented in the specialised scientific bibliography and published in journals and books specialised in the dentistry sector.

However, the most relevant are:

Post-operative discomfort – Localised pain – Localised inflammation – Haematomas – Haemorrhages – Localised or systemic infections – Difficulty speaking – Bone loss and/or fracture – Loss of the implant – Damage to adjacent teeth – Fractures to the implants and/or prosthetic components – Fracture of instruments – Injury to adjacent anatomical structures (inferior alveolar nerve or maxillary sinus).

Ziacom Medical SLU will update these side effects when necessary, based on the information obtained through clinical studies and post-sales information.

11. Information on responsibility, safety and guarantee

The indications for use and handling of the ZIACOM® products are based on the international published bibliography, the current clinical standards and our clinical experience with our products, and therefore should be understood as indicative general information. The handling and use of the ZIACOM® products, being out of the control of Ziacom Medical SLU, are sole responsibility of the user of them. Ziacom Medical SLU, its subsidiaries and/or official authorised distributors assume no responsibility, whether express or implicit, full or partial, for the possible harm or damages which may result from the improper handling of the product or from any other event not specified in its protocols and manuals for the correct use of its products.

The user of the product must ensure that the ZIACOM® product used is suitable for the intended procedure and purpose. Neither these instructions for use, nor the work and product handling protocols, absolve the user of this responsibility. The use, handling and clinical application of the ZIACOM® products must be carried out by qualified professional staff, with the necessary qualifications in compliance with the legislation in force in each country. The use, handling and/or application of the ZIACOM® products, whether fully or in part, at any stage, by unqualified staff or staff without the necessary qualifications, automatically renders any kind of guarantee void and can cause serious harm to the health of the patient.

The ZIACOM® products form part of their own system, with their own design characteristics and work procedures. This includes the dental implants, abutments or prosthesis components and the surgical or prosthetic instruments. The use of ZIACOM® products in

combination with elements or components from other manufacturers can cause failure of the treatment, cause tissue damage, cause damage to the bone structures, inadequate aesthetic results and serious harm to the patient's health. For this reason, only original ZIACOM® products must be used.

The clinician in charge of treatment is the sole responsible for ensuring that the original ZIACOM® products are used and that they are used in accordance with the corresponding instructions for use and working protocols throughout the orthodontic treatment process. The use of non-original ZIACOM® instruments used alone or in combination with any of the original ZIACOM® products will automatically void the guarantee of the original ZIACOM® products. Consult the Ziacom Medical SLU Guarantee Programme (available on the website or by contacting Ziacom Medical SLU, its subsidiaries or authorised distributors).

Warning. Not all the ZIACOM® products are available in all countries. Consult their availability.

ZIACOM®, Zinic®, Zinic®MT, Zinic®Shorty, OEX®, OEX®MT, OEX®ST, OIN®, OIN®MT, OIN®ST, OXP®, ZMK®, ZMR®, OST®, ZM1®, ZM1®MT, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, Galaxy®, ZV2®, Zinic®3D, TPlus®, XPlus®, ZPlus®, Z2Plus®, 3DPlus®, Kiran®, Kirator®, ZM-Equator®, Basic®, XDrive®, ZiaCam®, ZIACOR®, Tx30®, Zellplex®, ZellBone®, PlexGuide®, OsseosBCP®, OsseosTCP®, Osseolife®, Osseonova®, DSQ®, ZiTium® and all their logos are registered trademarks of Ziacom Medical SLU.

The ZIACOM® brand is a registered trademark, as are some of the products mentioned or not mentioned in this document. Ziacom Medical SLU reserves the right to modify, change, remove or develop any of the products, prices or technical specifications referenced in this document or in any of its catalogues without prior notice. All rights reserved. The whole or partial reproduction or publication of this documentation, by any means and in any format, is forbidden without the corresponding authorisation in writing from Ziacom Medical SLU.

12. Magnetic resonance

The safety and compatibility of the ZIACOM® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to see if they heat up or migrate in said environment.

13. Measurement function

The instruments with the function of measurement present pinpoint accuracy.

When using surgical drills, the length of the tip of the drill (+ 0.5 mm) must be taken into account as IT IS NOT INCLUDED in the measurements of the different laser brands, except in the Zinic®Shorty implants system.


14. Torque Wrench Maintenance

The torque wrench must be disassembled before cleaning by separating the individual parts. It is not necessary to use tools for this, only the adjusting nut must be completely unscrewed.

The loss of the plastic washer would negatively affect the accuracy of the instrument and should therefore only be removed in case of visible dirt. After cleaning, press the washer back in.

Lightly lubricate (before assembly) the indicated zones with white paraffinic oil without corrosion inhibitors or other additives (suitable for steam sterilisation and biocompatible). Use as little oil as possible in the marked areas.

We recommend an annual check of the torque wrench by the manufacturer.

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ZIACOM®, Zinic®, Zinic®MT, Zinic®Shorty, OEX®, OEX®MT, OEX®ST, OIN®, OIN®MT, OIN®ST, OXP®, ZMK®, ZMR®, OST®, ZM1®, ZM1®MT, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, Galaxy®, ZV2®, Zinic®3D, TPlus®, XPlus®, ZPlus®, Z2Plus®, 3DPlus®, Kيران®, Kirator®, ZM-Equator®, Basic®, XDrive®, ZiaCam®, ZIACOR®, Tx30®, Zellplex®, ZellBone®, PlexGuide®, OsseosBCP®, OsseosTCP®, Osseolife®, Osseonova®, DSQ®, ZiTium® und alle entsprechenden Logos sind eingetragene Marken von Ziacom Medical SLU.

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ZIACOM®, Zinic®, Zinic®MT, Zinic®Shorty, OEX®, OEX®MT, OEX®ST, OIN®, OIN®MT, OIN®ST, OXP®, ZMK®, ZMR®, OST®, ZM1®, ZM1®MT, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, Galaxy®, VZ2®, Zinic®3D, TPlus®, XPlus®, ZPlus®, Z2Plus®, 3DPlus®, Kiran®, Kirator®, ZM-Equator®, Basic®, XDrive®, ZiaCam®, ZIACOR®, Tx30®, Zellplex®, ZellBone®, PlexGuide®, OsseosBCP®, OsseosTCP®, Osseolife®, Osseonova®, DSQ®, ZiTium® ainsi que ses logotypes sont des marques déposées de Ziacom Medical SLU.

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ZIACOM®, Zinic®, Zinic®MT, Zinic®Shorty, OEX®, OEX®MT, OEX®ST, OIN®, OIN®MT, OIN®ST, OXP®, ZMK®, ZMR®, OST®, ZM1®, ZM1®MT, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, Galaxy®, VZ2®, Zinic®3D, TPlus®, XPlus®, ZPlus®, Z2Plus®, 3DPlus®, Kiran®, Kirator®, ZM-Equator®, Basic®, XDrive®, ZiaCam®, ZIACOR®, Tx30®, Zellplex®, ZellBone®, PlexGuide®, OsseosBCP®, OsseosTCP®, Osseolife®, Osseonova®, DSQ®, ZiTium® e todos os seus logótipos são marcas registadas da Ziacom Medical SLU.

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ZIACOM MEDICAL USA, LLC