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EN | INSTRUCTIONS FOR USE

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IMPLANT INSTRUCTIONS
FOR USE

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Ziacom®, Basic®, Di2gital Arch®, DSQ®, Galaxy®, Kiran®, Kirator®, OsseosBCP®, OsseosTCP®, PlexGuide®, Titansure®, Titansure® Active, Tx30®, XDrive®, Zellplex®, ZiaCam®, Ziacom®3D, Ziacor®, Zinic®, Zinic®MT, Zinic®Shorty, ZiTium®, ZM1®, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, ZM-Equator®, ZMK®, ZMR®, ZPlus®, ZV2® y todos sus logotipos son marcas registradas de Ziacom Medical SL.

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Ziacom® orthodontic implants

INSTRUCTIONS FOR USE

Instructions for use applicable to the DSQ® orthodontic implants systems.

IMPORTANT INFORMATION

Read these instructions for use carefully before using the product.

This document contains basic information for the use of original Ziacom® orthodontic implants, micro-implants or mini-screws, hereinafter “Ziacom® orthodontic implants” or simply “Ziacom® products”. This document has been written as a quick reference guide for the physician in charge of the treatment of orthodontic, henceforth, “User”. It does not provide sufficient indications and technical specifications for the use of Ziacom® products. It is neither an alternative nor a substitute for specialised training and professional clinical experience.

Ziacom® products must be used with proper treatment planning and following the working protocols established by Ziacom Medical SL. Read carefully: the specific working protocols for this product, as well as the instructions for use and the instructions for cleaning, disinfection and sterilisation before using this Ziacom® product. You can consult them on our website www.ziacom.com or request them from your official distributor authorised by Ziacom Medical SL.

1. Description of the system

Ziacom® products consist of orthodontic implants, dental implants, prosthetic abutments or components and surgical and prosthetic instruments.

Ziacom Medical SL has created two different types of orthodontic implants to serve as temporary anchor points in orthodontic treatment according to the existing scientific literature and current clinical standards. Ziacom Medical SL develops, manufactures and commercialises latest-generation implant products and a wide range of restorative solutions to respond to the particular needs of each patient.

2. Description of the products

Orthodontic implants

Ziacom® orthodontic implants are manufactured in biocompatible grade 5 titanium alloy, in accordance with ASTM F136-13 and UNE-EN-ISO 5832-3:2012 standards, without surface treatment to avoid excessive osseointegration that hinders explantation at the end of treatment. Ziacom® orthodontic implants are manufactured in different diameters, lengths, heads (cross or button type) and collars depending on their intended purpose or location.

3. Indications for use

It is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and study of models, as well as to determine the possible contraindications for the treatment, prior to carrying out a treatment with Ziacom® products. A correct orthodontic treatment that requires the use of orthodontic implants must start from a diagnosis, planning and design of orthodontic treatment that meets the needs and expectations of the patient.

Orthodontic implants

Ziacom® orthodontic implants are indicated as a temporary skeletal anchorage point for orthodontic treatment in both the jaw and maxilla, and are intended to be attached to different orthodontic appliances to generate or support desired tooth movements. They are supplied with different head designs and in different lengths and diameters so that they can be adapted to different clinical cases.

Ziacom® implants are single-use products. The head design of orthodontic implants allows their use together with other products and accessories commonly used in orthodontics (orthodontic wires, steel ties, elastic ligatures, brace chains, springs, etc.), for more information consult the catalogue.

The insertion and removal of the orthodontic implant must be performed by a qualified professional using the recommended instruments and following the indications in the catalogues and working protocols.

Important

Carefully consult the clinical indications for each orthodontic implant model in the relevant catalogues and working protocols before planning treatment. Do not use Ziacom® orthodontic implants in locations or for purposes that are not described as indications in your catalogues or working protocols.

4. Contraindications

A preoperative medical examination of the patient is necessary to determine whether there are absolute and/or relative contraindications or risk factors for treatment with orthodontic implants. It is the user's responsibility to evaluate the potential benefits and risks of using orthodontic implants in orthodontic treatment and to make appropriate decisions on a case-by-case basis. Ziacom® products should not be used in patients who lack the minimum medical conditions for orthodontic treatment and restoration with orthodontic implants.

General contraindications:

Patients with systemic diseases (disorders in blood coagulation, uncontrolled endocrine disorders, rheumatic diseases, liver diseases, neoplasms, diabetes or other serious diseases), psychological or psychic alterations, poor oral hygiene, severe or uncontrolled periodontal disease, inadequate quantity/quality of bone that prevents the retention of the orthodontic implant in the bone cortex, inadequate oral habits (aggressive tooth brushing, smoking, foods with high sugar content, oral piercings...).

Ziacom Medical SL will update these contraindications when necessary, based on the information obtained through clinical studies and post-sales information.

5. Storage

Ziacom® products must be stored in their original packaging, in a clean, dry place, protected from sunlight and adverse conditions. The storage temperature must be lower than +55°C and higher than -10°C.

6. Packaging and supply state

Ziacom® orthodontic implants are supplied in a cardboard box with a product identification label, which acts as a security seal against handling.

This package contains:

- This explanatory product information leaflet in paper or electronic format.
- Thermo-sealed plastic blister with removable labels on the back, indicating the manufacturer, complete product reference, lot number, dimensions of the orthodontic implant and the relevant symbol (explained in these Instructions for use). Inside is the orthodontic implant.

The removable labels serve as a traceability control of the implant. This is achieved by sticking one label to the implant card for the patient, and the rest for the clinic's record.

Ziacom Medical SL guarantees that all its products follow an exhaustive process of manufacture, control and cleaning before being packaged. Before using Ziacom® products you must check the correct integrity of the package and make sure that it is not damaged. If there is any defect or damage, it must not be used and must be notified to Ziacom Medical SL or to the official authorised distributor.

7. Cleaning and sterilisation

Ziacom® orthodontic implants are supplied non-sterile.

Orthodontic implants that will come into contact with the mouth must be cleaned, disinfected and sterilised prior to use following the cleaning, disinfection and sterilisation protocol for orthodontic implants, instruments, surgical boxes and abutments of Ziacom Medical SL. You can consult it in our web page www.ziacom.com or in the catalogues of Ziacom Medical SL.

8. Precautions

Insertion and instruments

The working protocols describe in detail the important precautions and factors to be taken into account during the orthodontic implant placement process. Before insertion, it is essential to consult the working protocol of the product, as well as all indications for its correct use in the Ziacom® product catalogues.

Special attention must be paid and the necessary measures must be taken to prevent the orthodontic implant or instruments from being swallowed or bronchoaspirated by the patient.

Pay special attention not to exert excessive forces or torques during the insertion of the orthodontic implant. If during insertion you detect a variation in bone consistency, stop insertion immediately and make sure that you are not inserting into inappropriate bone or dental structures such as dental roots.

Attention

The maximum insertion depth of orthodontic implants is up to the collar of the orthodontic implant. Avoid applying excessive insertion torque to orthodontic implants as this can cause serious damage to the orthodontic implant and bone. The application of the insertion torque should be under the responsibility of the user.

9. Warnings

Treatment planning and placement of Ziacom® orthodontic implants require specific dental training. Product specifications alone do not guarantee proper use. Users are encouraged to take specialised theoretical and practical training courses to understand the appropriate working techniques and protocols for each product, including biomechanical, radiographic and prosthetic requirements associated with orthodontic treatments.

It is the user's obligation to know and be informed of the current state of the technique for each product and its possible applications.

Before using Ziacom® products it is necessary to know and be familiar with the corresponding working procedures. The patient must have sufficient anatomical and psychological conditions to undergo dental treatment. The user must ensure the correct planning of the treatment, guaranteeing its safety with a minimum margin of error and respecting the vital structures of the mouth and the general health of the patient.

The procedure for the use of the insertion instruments required for the placement of orthodontic implants is specified in the working protocols for each type of orthodontic implant. Orthodontic implant placement and orthodontic planning should be tailored to individual patient conditions. The point of insertion should be chosen so that the load should have a 90° angle or be perpendicular to the longitudinal axis of the orthodontic implant. Under or overloaded treatment, an inappropriate choice of size or an inappropriate position to support the expected loads can lead to failure of the orthodontic implant or treatment (overload or fatigue).

A qualitatively or quantitatively insufficient alveolar bone, the appearance of infections or diseases in general, and changes in the habits of the patient, are some potential causes of failure of the treatment. Lack of bone or soft tissue can lead to poor outcome and inadequate loading can lead to failure of the entire restorative treatment.

The reuse of single-use products carries the danger of possible deterioration of their properties. The geometry of the product may have suffered damage during previous use or due to an unsuitable use. This implies risk of orthodontic failure and/or damage to the patient's health.

Electrosurgery is not indicated due to the conductivity of orthodontic implants.

10. Adverse or side effects

Many adverse or undesired effects can present themselves in orthodontic. These are documented in the specialised scientific bibliography and published in journals and books specialised in the dentistry sector.

However, the most relevant are:

Operative or postoperative pain, inflammation/infection of the mucosa of the area, haematomas or haemorrhages, bone mobility, loss or fall of the

orthodontic implant, alveolar fracture, fracture of the orthodontic implant, damage to adjacent roots or teeth, injury to adjacent anatomical structures (nostrils, maxillary sinus, vascular/nerve structures).

Ziacom Medical SL will update these side effects when necessary, based on the information obtained through clinical studies and post-sales information.

11. Information about responsibility, safety and guarantee

The handling and use indication of Ziacom® products are based on published international literature, current clinical standards and our clinical experience with our products and should therefore be understood as indicative general information. The handling and use of Ziacom® products, as they are beyond the control of Ziacom Medical SL, are the sole responsibility of the user. Ziacom Medical SL, its subsidiaries and/or its authorised official distributors decline all responsibility, express or implicit, total or partial, for the possible damages caused by the bad handling of the product or by any other fact not contemplated in its protocols and instructions for the correct use of its products.

The user of the product must ensure that the Ziacom® product used is suitable for the procedure and the intended purpose. Neither these instructions for use, nor the working protocols of the products, exempt the user from this obligation. The use, handling and clinical application of Ziacom® products must be carried out by qualified professional personnel with the necessary certification according to the legislation in force in each country. The use, manipulation and/or application, in whole or in part, of Ziacom® products in any of their manufacturing phases by unqualified personnel or without the necessary qualifications, automatically voids any type of guarantee and may cause severe damage to the health of the patient.

Ziacom® products have their own design features and working protocols, so the insertion of orthodontic implants must be done with the adapters and insertion instruments recommended by Ziacom Medical SL in the catalogue and in the working protocols. However, the design of orthodontic implants enables them to be used together with other products and accessories commonly used in orthodontics (orthodontic wires, steel ties, elastic ligatures, brace chains, springs, etc.). Insertion of Ziacom® orthodontic implants with adapters or instruments from other manufacturers may result in insertion failure, deterioration of the orthodontic implant, injury to the insertion site by causing tissue damage, bone structure damage, inadequate aesthetic results and serious damage to the patient's health. For this reason, insertion should only be performed with instruments and adapters recommended by Ziacom Medical SL.

The clinician in charge of treatment is the sole responsible for ensuring that the original Ziacom® products are used and that they are used in accordance with the corresponding instructions for use and working protocols throughout the orthodontic treatment process. The use of non-original Ziacom® instruments used alone or in combination with any of the original Ziacom® products will automatically void the guarantee of the original Ziacom® products. Consult the Ziacom Medical SL Guarantee Programme (available on the website or by contacting Ziacom Medical SL, its subsidiaries or authorised distributors).

Warning

Non entire Ziacom® products are available in all countries. Please consult for availability.

Ziacom®, Basic®, Di2gital Arch®, DSQ®, Galaxy®, Kiran®, Kirator®, OsseosBCP®, OsseosTCP®, PlexGuide®, Titansure®, Titansure® Active, Tx30®, XDrive®, Zellplex®, ZiaCam®, Ziacom®3D, Ziacor®, Zinic®, Zinic®MT, Zinic®Shorty, ZiTium®, ZM1®, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, ZM-Equator®, ZMK®, ZMR®, ZPlus®, ZV2® and all their logos are registered trademarks of Ziacom Medical SL.

The brand Ziacom® is a registered trademark, as some of the products mentioned or not in this document. Ziacom Medical SL reserves the right to modify, change, eliminate or evolve any of the products, prices or technical specifications referenced in this document or any of its catalogues without prior notice. All rights reserved. Reproduction or publication in whole or part of this documentation, in any medium or format, without the corresponding written authorization of Ziacom Medical SL is prohibited.

12. Magnetic resonance

The safety and compatibility of Ziacom® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to see if they heat up or migrate in said environment.

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Ziacom®, Basic®, Di2gital Arch®, DSO®, Galaxy®, Kiran®, Kirator®, OsseosBCP®, OsseosTCP®, PlexGuide®, Titansure®, Titansure® Active, Tx30®, XDrive®, Zellplex®, ZiaCam®, Ziacom®3D, Ziacor®, Zinic®, Zinic®MT, Zinic®Shorty, ZiTium®, ZM1®, ZM4®, ZM4®MT, ZM8®, ZM8®N, ZM8®S, ZM-Equator®, ZMK®, ZMR®, ZPlus®, ZV2® e todos os seus logótipos são marcas registadas da Ziacom Medical SL.

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Symbols of labelling codes and used in catalogues, standard ISO 15223-1

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MDD CE certification and notified body

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Caution, Consult accompanying documents

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Rx Only

Caution: federal law prohibite dispensir without prescription.

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
EN Once the instruments have been sterilised, use the insertion key in the handle or the counter-angle insertion key to hold the implant for clinical use.

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EN Remember to remove the identification labels from the implant to stick them on the Implant Card and on the patient history sheet in order to document the traceability of the reference and batch of the implant.

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